

MediaLINC is a powerful software platform providing patients with access to educational materials and entertainment options through their in-room television. It is an integrated solution for the delivery of digital patient engagement services such as information, education, entertainment and internet to the patient bedside using a set top box and standard HD TV



MediaLINC features a unique user interface designed for the hospital environment. It can be used in single occupancy rooms as well as multi-bedded rooms or public areas. Utilizing the power of the LINC OS platform it can deliver a range of services:

MyCare

MyCare consists of a range of patient focused applications designed to improve the hospital experience:

Hospital Information

An easily customer maintainable range of information useful to the patient.

Meal Ordering

An intuitive, customizable interface that facilitates food ordering for patients and visitors.

Surveys

Easy-to-use intuitive interface for patient satisfaction surveys and educational quizzes.



Entertainment

MediaLINC features a range of entertainment options available via simple to use intuitive interfaces designed to match the services being used.

All services can be controlled, managed and, where required, packaged by the customer using the powerful LINC OS platform.



TV – Live TV Streaming from a range of channels.

Radio – Live streaming of online radio stations

Video on Demand – Access to a library of content from a range of content providers

Talking Books – a library of audio books

Music Library – A centralized Jukebox


Note: All content remains in the control of the customer and may be subject to additional subscription costs

Applications

MediaLINC allows any TV compatible Android application to be deployed on the solution.

Internet Access

MediaLINC provides internet access via an HTML 5 compatible browser supporting direct entry of URLs as well as easily customizable shortcuts. For casual use, an onscreen virtual keyboard is available, but for ease of use we can supply a wireless keyboard with trackpad.



Settings

There are a variety of settings that help users get the most out of MediaLINC. These include:

Language Support – The patient can change the language of the user interface.

Automated Background – The interface automatically switches to a softer, darker background at night to be more soothing.

Simplified View – There is an option to simplify the interface by removing items such as animation and background, that may be confusing for some patients.

Media Share

If used in conjunction with PatientLINC, the user can show what is being viewed on the bedside smart device on the TV.

Billing

MediaLINC supports pre-paid scratch card billing to allow the hospital to create a valuable new revenue stream. The hospital maintains complete control over pricing and bundling. Additional billing methods will be introduced in future releases.



Technical Specifications

Set Top Box

The MediaLINC solution uses a lightweight android-based set top box that automatically ensures you are running the latest version of the software through over-the-air updating.

Specifications: Interfaces: HDMI, 3xUSB, Optical, LAN, AV, TF Card Slot

Power Supply: 5V/2A

Dimension: 150x115x27mm

Net Weight: 0.2 Kg

Connectivity Requirements

MediaLINC requires fixed CAT5/CAT6 network to each location.

TV Compatibility

MediaLINC has been designed to work with any HD TV with HDMI connectivity. Owing to differences in TV models, Lincor reserve the right to test against selected TVs prior to final deployment.